



The MR Conditional folding portering chair is designed to be used in an MR environment and is manufactured & tested upto 3 Testa

Wheelchair for your MRI Facility



Produced from non-ferromagnetic quality materials with high quality and patient comfort. The product is non affected by magnetic field.



The MR Wheelchair can carry weight upto 120 kg and is specially designed for use inside 3T MRI rooms.

3T



The MR Wheelchair is made for easy handling for movement anywhere inside MRI room



Kryptonite Solutions,
Plt.32, Behind Suvarna Hospital,
Shimpoli Road, Borivali West,
Mumbai - 4 00 092, India
Phone: +91 22 28995246
Email: info@kryptonite.global